

Chromascan Modules	Application Dyes
Chromascan UV365 pack module EX365nm-EM590nm	Qdot 565, Qdot 655, Qdot 705, TLC plates, microplate, Alexa 350, DAPI
Chromascan Light Blue pack module EX440nm-EM500nm	DAPI, CFP, Cerulean, Alexa Fluor 405, Cascade Blue, Pacific Blue, DyLight 405, Atto 425
Chromascan Blue pack module EX480nm-EM550nm	FITC, Alexa Fluor 488, GFP, SYTOX Green, Fluorescein, Cy2, Sypro Ruby, DyLight488
Chromascan Deep Blue pack module EX480nm-EM600nm	YFP, eYFP, Venus, Alexa Fluor 514, FITC, mCitrine
Chromascan Green pack module EX540nm-EM600nm	Rhodamine, Alexa Fluor 532, Alexa Fluor 555, Cy3, PE, TRITC, ProQdiamond, DyLight 549, mRuby3
Chromascan Deep Green pack module EX540nm-EM650nm	Cy3.5, Atto 565, Rhodamine 6G
Chromascan Orange pack module EX580nm-EM650nm	DsRed, mCherry, Cy3.5, Alexa Fluor 568, Texas Red, Atto565, Atto594, Alexa594, mStrawberry, mKate2
Chromascan Red pack module EX640nm-EM700nm	Alexa Fluor 647, Alexa Fluor 660, Cy5, APC, Atto633, Atto 647N, DyLight 633, DyLight 650
Chromascan NIR pack module EX680nm-EM750nm	Alexa Fluor 680, Alexa Fluor 700, Cy5.5, IRDye 680RD, Atto 680, Atto 700, APC-Cy7, DyLight 680
Chromascan LIGHT IR pack module EX740nm-EM750nm	Alexa Fluor 750, Cy7, IRDye 800CW, VivoTag-S 750, DyLight 755
Chromascan IR pack module EX740nm-EM800nm	Cy7, IRDye 800RS, ZW800-1, Atto 740, HiLyte Fluor750
Chromascan FAR IR pack module EX780nm-EM850nm	Alexa Fluor 790, Cy7.5, IRDye 800CW, VivoTag-S 800, DyLight 800